Dräger Alcotest[®] 6510

Atemalkohol-Messgerät
Breath Alcohol Monitor

de

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Seite 2

en

Instructions for Use

Page 25

Inhalt

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Content

For Your Safety	26	Quick Menu	35
Intended Use	26	Main Menu	35
What is What?	27	Error	41
Case	27	Maintenance	41
Measuring instrument	27	Maintenance intervals	41
Preparation	28	Storage	41
Inserting the Batteries	28	Cleaning	42
Attaching the brief instructions for use	29	Changing Batteries	42
Fixing the hand strap	29	Calibration	43
Symbols and readings	29	Disposal of the instrument	43
Fitting the Mouthpiece	30	Fault, Cause, Remedy	44
Switching ON and OFF	30	Technical Data	45
Performing the Measurement	31	Order List	47
Conditions	31		
Requirements of the person to be tested	31		
Automatic Measurement	31		
Errors in delivering the breath sample	32		
Removing the mouthpiece	33		
Additional possibility for sampling	33		
Passive sampling	33		
Manual sampling	33		
Menu	34		
Activation	34		
Operation	34		
Symbols on the Display	34		

For Your Safety

Strictly Follow the Instructions for Use

Any use of the instrument requires full understanding and strict observation of these instructions. The instrument may only be used for the purposes specified here.

Maintenance

The instrument must be inspected and serviced on a regular basis by qualified personnel. Repair and general overhaul of the instrument may only be carried out by trained service personnel.

We recommend signing a service contract to have Dräger carry out all maintenance jobs.

Only Dräger original parts may be used for maintenance.

Be sure to read the information contained in the chapter "Maintenance" on page 41.

Accessories

Only use accessories mentioned in the order list.

Safety Symbols Used in these Instructions for Use

While reading this manual, you will come across a number of warnings concerning some of the risks and dangers you may face while using the instrument. These warnings have signal words to make you aware of the expected degree of danger. These words, and the hazard they describe, are as follows:

WARNING

Indicates a potentially hazardous situation which, if no safety precautions are taken, could result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation which, if not avoided, could result in physical injury to you, or damage to the product. It may also be used to alert against unsafe practices.

NOTICE

Additional information on how to use the instrument.

Intended Use

For quick measurement of a person's alcohol concentration through measurement of the breath alcohol concentration. For mobile use with batteries.


Preparation

Dräger Alcotest® 6510 is supplied with standard calibration by the manufacturer. If the calibration date exceeds the period of 6 months: recalibrate instrument or have it calibrated by DrägerService. The date of the last calibration can be found in the quick menu (see page 35).

® Alcotest is a registered trademark of Dräger.


What is What?

Case

- 1 Dräger Alcotest® 6510 monitor
- 2 Mouthpieces, 3 pieces
- 3 Batteries, 2 pieces
- 4 Hand strap
- 5 Transport case
- 6 Brief instructions for use

Measuring instrument

- 1 Display
- 2 Menu button "UP"
- 3 "OK" button
- 4 Menu button "DOWN/MENU"
- 5 Socket
- 6 Lamp red/green
- 7 Mouthpiece receptacle

Preparation



Inserting the Batteries

CAUTION


Charging batteries while in the instrument will cause damage to the instrument.

Upon initial operation or if the battery symbol on the display indicates **"empty"**.

- Gently press the lock on the battery cover and remove the cover.
- Observe correct polarity when inserting the batteries.
- Close the battery cover.


Batteries to be used
2 x 1.5 V alkaline cells
(Mignon, LR6, AA)
or
rechargeable battery
2*1.2 V NiMH

Attaching the brief instructions for use


An area for the "brief instructions for use" sticker is provided on the back of the instrument.

- The instrument is delivered with the brief instructions in English, German, French and Spanish language.
- Remove sticker from the carrier foil.
- Position the sticker onto the area provided on the back of the housing.

Fixing the hand strap

- Insert the smaller loop into the hole in the housing.
- Thread the longer loop through the eyelet and pull through.

Symbols and readings

The measured breath alcohol value is given in mg/L for the basic setting of the instrument. The instrument can be switched over to all usual measuring units (see page 38).

Symbols

In these instructions for use several symbols are used corresponding to the messages indicated on the display:

Battery Symbol

Indicates the state of charge of the batteries:

Full battery


50 % battery capacity

Low battery

NOTICE

The symbols under the battery flap mean: Attention! Please consult the user instructions!
Before changing batteries, please read the instructions. The service socket is not a charging socket. Never insert a power cable!

Buttons

"OK" button

Confirms entries or selected functions.
Button to switch instrument on and off.


"DOWN" and "UP"


Setting of values and selection of menu functions.

Fitting the Mouthpiece


- Open packaging. To keep the mouthpiece hygienic, retain the portion to come in contact with mouth covered by packaging until mouthpiece is correctly fitted to the instrument.
- With the grooved markings facing up, press the mouthpiece into mouthpiece receptacle until audible sound confirms correct engaging.
- Remove packaging completely and commence breath test.
- The mouthpiece can be used with left or right orientation.

0323757_1_en.apa

0323757_1_en.apa

0323757_1_en.apa

Switching ON and OFF

Switching ON

- Keep  - button pressed for approx. 1 second until start window is displayed.

01323757_1_en.apa

Switching OFF

- Keep  - button pressed.
- The instrument switches off after 2 seconds.

01423757_1_en.apa

Automatic Switch OFF

The instrument switches off automatically 4 minutes after the last use.

When the batteries are discharged the instrument switches off automatically.

0323757_1_en.apa

Performing the Measurement

Conditions

- Use a new mouthpiece for each test subject.
- Do not apply too high alcohol concentrations on to sensor, e.g. by rinsing the mouth with high-grade alcohols shortly before the measurement. By doing so the sensor service life is reduced.
- Keep sufficient distance from transmitters and mobile phone antennas.

Requirements of the person to be tested

- Maintain a waiting period of at least 15 minutes after drinking alcohol!
Residual alcohol in the mouth can distort the measurement. Aromatic drinks (e.g. fruit juice), mouth sprays containing alcohol, medicines can interfere with measurements. Rinsing the mouth with water or non-alcoholic drinks does not reduce this waiting period!
- Before sampling the person to be tested must breathe normally. Quick inhaling and exhaling through the mouth must be avoided.
- The person to be tested must be able to provide the required minimum air volume of 1.2 L (factory setting). The breath flow must be constant for a certain minimum blowing period (company setting 2 seconds). Volume and time can be changed via the menu, if necessary (see page 37).

Automatic Measurement

The sampling is activated automatically after the minimum breathing volume and the minimum blowing period are reached.

Preparation

- Insert new mouthpiece in the holder.
(see page 30).
- Switch on the instrument (see page 30).

After 2 seconds, the display indicates:

- the user info screen
(if activated, see page 37)
or
- **"WAIT"**.
- The user info screen can also be viewed by pressing the - button. Close with the - button

After approximately 6 seconds **"READY"** is displayed and a short audible signal is given. Simultaneously, the current test number is displayed. The green lamp flashes.


The instrument is ready to start the measurement.

Performing the Measurement

- Ask the test subject to blow evenly into the mouthpiece without interruption. A continuous tone sounds and the green lamp flashes when a sufficient breathing flow volume has been reached.

While the sample is taken the display shows:
"Blow"


When the sampling was successful the green light will cease and the continuous tone will stop.

The display shows:
"WAIT Analyzing"

Test Result

After 5 to 25 seconds (depending on the instrument temperature and the measured concentration), the measuring result appears on the display.

The value is shown with the respective measuring unit.

The standard setting by the manufacturer is mg/L (milligram ethanol per litre breath air).

The setting can be changed via the main menu (see page 38).

Errors in delivering the breath sample

If the sampling volume is too low the display shows: **"Insufficient Volume"**.

The horn sounds briefly and a red light flashes.

To repeat press the  - button.

The instrument can be used again after approx. four seconds.

When exhaling improperly, e.g. by abruptly stopping to exhale or sucking at the end of the sampling, the display shows:

"Blow Interruption"

A horn sounds briefly and the red lamp flashes.


To repeat press the  - button. Time period until the instrument is ready for measurement again: 4 to 30 seconds.

Further Measurement

- By pressing the  - button a new measurement is prepared.

During the recovery period of the sensor (see "Technical Data" on page 45), the display shows: **"WAIT"**

NOTICE

Temporarily switching off the instrument will not reduce the waiting period!

Removing the mouthpiece

- Pull the mouthpiece out of the mouthpiece holder.
- Dispose of the mouthpiece according to the local regulations.

Screen 1:
No alcohol available
or:

**NO
ALCOHOL**

02523757 1 en.ind

Screen 2:
In case of alcohol
● Press the  - button for a new measurement.

ALCOHOL

02523757 1 en.ind

Additional possibility for sampling

The functions can be switched on or off in combination also in the main menu (see page 35).

Passive sampling

During the passive sampling the sampling system is manually activated by the user of the instrument. Possible applications for this function are testing of the ambient air or the exhaled air for alcohol concentrations of persons not able to cooperate (e.g. unconscious persons).

- Hold the **instrument without a mouthpiece** into the air to be analysed.
- Press the  - button shortly.

Manual sampling

The sampling can also be activated manually, if the person to be tested is not able to provide the required minimum volume.

Preparation as for the automatic measurement (see page 31).

- If "Blow" is indicated on the display after exhaling has started, briefly press the  - button shortly.

When the measured value is displayed
"Manual" is indicated additionally.

0.00
Manual mg/L

02523757 1 en.ind

Menu

Activation

Turn the instrument on. When **"WAIT"** or **"READY"** appears on the display the menu functions can be called up.

- For quick menu press button  shortly.
- For main menu keep button  pressed for at least 3 seconds.

Operation

- Select functions using  and  buttons.
- The selected function is highlighted.
- Activate the selected function by pressing  - button.

When the selected function is displayed the top line contains the current menu function.

- Selection " **return to measurement**" or " **return to menu**" and pressing the  - button brings the instrument into the measuring mode or to the former menu level back.

If no button is pressed within 120 seconds, the instrument automatically returns to the measuring mode.

Symbols on the Display

Menu

-  Return to next higher menu level
-  Menu item with sub-menu
-  Open menu item
-  Selected function
-  Further items below the visible items
-  Further items above the visible items


Function

-  Activated item
- 1/6** Actual screen (1) of all screens (6)

Quick Menu

■ Last 10 tests

Selection of the tests with ▼ and ▲ .
Back with (OK) - button.

■ Last cal. date

Display of last calibration date. Back with (OK) - button.

■ Select language

English, German, French, Spanish, user defined language to be load in the instrument.


Set contrast

Change by pressing ▼ and ▲ . Back with (OK) - button.

Main Menu

The main menu consists of three sub-menus

- Information
- Settings
- Calibration


Access to the main menu is protected by a PIN.

Entering the PIN

The display shows:


"Enter PIN".

- Press (OK) - button.

- Select number or letter by pressing ▼ and ▲ .

- Acknowledge by pressing (OK) - button and enter next digit.


After the fourth entry the main menu is displayed.

It can be changed under settings **"Change PIN"** (see page 37).

In case of a wrong PIN the display shows **"Wrong PIN"** together with the entered PIN.


- (OK) Press - button and enter new PIN.

Information

The menu item "**Information**" includes the following items:

- Last 10 tests
- Last cal. date
- Configuration
- Test number
- Software version

Additional sample

Activated additional functions such as manual or passive sampling.

Volume

Minimum exhalation volume in litres.

Time

Minimum blowing time in seconds.

Test Number

Display of the number of the next test.

Software Version

Display of the part number and the software version.

Last 10 tests, last cal. date

see quick menu (see page 35).

Configuration

Under configuration you will find the current instrument settings.


These settings can be changed under the menu item settings (see page 37). Here you will find more information concerning these items.

Display mode

Display of the measured value in numbers or as a range.

Measurement Unit

- Measuring unit of the alcohol concentration (factory setting mg/L).

Pass, Alert, Fail, High

Set limit values for the display modes PAF, PF and ZPWF.


Auto Sample

Type of sampling: in case of definite volume or end of breath sample.

Settings

The menu item **"Settings"** includes the following items:


- Instrument
- Measurement

+ Instrument

■ Select language, set contrast

see quick menu (page 35).


■ Test number

Activate or deactivate the display of the test number.

■ Change PIN

Enter the new 4 digit PIN by pressing ▼ , ▲ and OK - buttons.

(Numbers and letters are possible).

NOTICE

If a changed PIN code is no longer available, only DrägerService has access to the menu.



■ Individual information

Entry of up to 3 lines with 15 characters each pressing ▼ , ▲ and OK - button (Numbers and letters are possible).

For finishing a line keep OK - button pressed until the line appears inverted.



Confirmation with OK - button stores the text.

These symbols are indicated when switching on the instrument.

+ Measurement

Setting of parameters and limit values for the measurement.

Display mode


Selection of the result presentation:

a) numerical = corresponding to the set measuring unit

b) Allocation of the result to in concentration ranges with selected range limits.

Available:

- PAF = Pass, Alert, Fail
- PF = Pass, Fail
- ZWPF = Zero, Warn, Pass, Fail

Measuring Unit

Selection of the measuring unit for the final result:

%, ‰, mg/L, µg/L, g/210 L, µg/100mL, µg/%, g/L, mg/100mL, mg%, mg/mL

Result resolution

Selection of the measurement resolution.

Sample Mode

Auto sample

Selection of the mode of time for sampling:


a) definite volume

Sampling is complete when minimum volume is reached.

or

b) end of breathing


Sampling is complete at the end of exhaling.

Additional sample

Manual activation of sampling and lower minimum volume

- Passive
- Manual


Also refer to page 33 and page 33.

BAC - Limits

Three range limits "Pass, Alert, Fail" for the display modes PAF, PF and ZWPF.

Select and enter by pressing  ,  and  --button.

The illustration on the following page shows the connection between the limits and the display of the result.

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
Blow criteria

Setting of the minimal breathing volume and the minimum blowing time for a valid sample. After selection of the function, entry in litres (L) by pressing ∇ , \triangle and

Factory setting is 1.2 L and 2 seconds.

Control Mode

Increased resolution with 3 decimal places only for measuring in mg/L and µg/L.


Parameter settings are reset after switching off the instrument.

It is possible to activate all functions simultaneously.

Calibration

■ Last cal. Date

In the sub-menu item "**Last cal. date**" the calibration date last entered can be viewed,
see quick menu (see page 35).


+ Options

Setting of the following items with the  ,
 and the  - button.

■ Setpoint

concentration of test gas in the set measuring unit

- calibration
- accuracy check

(permissible calibration range, see "Technical data" page 46).

■ Gastype


wet = bath gas, dry = dry gas

■ Altitude

height above sea level,
only necessary for dry test gas
(selectable unit: metres or feet)

■ Calibration

- After pressing the  -button the display shows the "**Setpoint**" and the "**Gas type**".
- After pressing the  -button the instrument will be ready for calibration.
- When "**Supply gas**" is shown supply respective gas.
- A minimum flow rate of 3 L/min to obtain "**Gas flow OK**" on display.
- The instrument performs the calibration automatically.
- After successful calibration the message "**Calibration OK**" is displayed.
- A message "**Calibration Error**" displayed requires the calibration to be repeated.
Repeat calibration in case.




Accuracy check

- Enter the **"Accuracy check"** item by using the **(OK)** -button. The display shows the calibration gas to be used.
- Follow the sampling procedure described for calibration.

After successful accuracy check the result is displayed.

Enter calibration date


- Use **(↓)** , **(↑)** and **(OK)** - button in **"Enter cal. date"** menu to enter the format and the date.

Error

An error message indicates malfunctions of the instrument.

The number below the message is for service purposes. If the failure still occurs after the instrument is switched on and off several times, contact DrägerService.

Maintenance

Maintenance intervals

As required:

- Changing batteries, see page 42.

Every six months:

- Check calibration, see page 43.

Storage

- The instrument consumes a small amount of current even when switched off. Therefore, remove the batteries if the instrument will not be used for a longer period of time.

NOTICE

If case of long-term storage

Remove the batteries from the instrument if it will not be used for more than six months. This prevents possible damage to the instrument caused by leaking battery acid.

⚠ WARNING

Danger of explosion!

Do not throw used battery into open fire and do not try to force them open.


Disposal of batteries in accordance with national regulations.

Cleaning

- Do not submerge in liquid
– liquid must not enter the sockets!

To clean:


- wipe with a disposable tissue.
- do not use cleaning agents containing alcohol or solvents.

Changing Batteries

If the battery symbol on the display indicates "empty", change the batteries.

- Gently press the lock on the battery cover and remove the cover.
- Remove used batteries and insert new ones – consider correct polarity.
- Close the battery cover.

Required batteries:

2 x 1.5 V alkaline cells (Mignon, LR6, AA)

or

rechargeable battery 2*1.2 V NiMH

- For continuity of use it is recommended always to have replacement batteries or an additional battery supply available.

⚠ WARNING

Danger of explosion!

Do not throw used battery into open fire and do not try to force them open.

Disposal of batteries in accordance with national regulations.

Calibration

- Only by authorised personnel.
- Calibration recommended every 6 months.
- in sub-menu item **"Last cal. date"** the calibration date last entered can be viewed.

Preparation measures

- Call up **"Calibration"** in the main menu and select menu item **"Options"**.
- Set the **"Setpoint"**, the **"Gas type"** and the **"Altitude"** in menu item **"Options"**.

NOTICE

The settings of gas type and height above the sea level are important only for the calibration.

Calibration using the Dräger-Mark II Simulator

Follow instructions for use of the simulator.

Prepare simulator for operation.

Connect Dräger Alcotest® 6510 mouthpiece via adapter with the simulator outlet. To avoid condensation the connection must not be longer than 5 cm.

Calibration with Dry Gas

Connect instrument with test gas bottle via a hose (maximum length 1 m).

Performing the Measurement

Procedure as calibration described in the menu on page 40.

- After successful calibration an **"Accuracy check"** for testing the instrument must be performed.
 - The settings gas type and altitude remains.
- The result of the accuracy check shall be compared with the set point.

Disposal of the instrument

Since August 2005, regulations concerning the disposal of electrical and electronic instruments have come into force in the entire European Union; these regulations are specified in the Directive 2002/96/EC of the European Parliament and in national laws relating to this instrument.

Special collection and recycling possibilities have been established for private households. However, since the instrument has not been registered for use in private households, it cannot be disposed of as special household waste. Therefore, you should send it to your national Dräger distribution organisation for disposal, which you can contact for further information.

Fault, Cause, Remedy

Fault	Cause	Remedy
Instrument switches off automatically.	batteries discharged.	Replace the batteries, page 42.
	The instrument has been ready to measure for more than four minutes.	
Instrument leaves menu-mode automatically and shows measuring mode.	Instrument has an automatically switch off of the menu after 120 sec. without use.	Reselect menu item.
Instrument does not power up.	Batteries are discharged.	Insert new batteries.
Insufficient volume/sample invalid.	The test subject is not breathing continuously or vigourously enough.	Ask to blow harder and more evenly.

Fault in display	Cause	Remedy
"Calibration error"	Calibration not successful.	Repeat accuracy check, if fault reappears, call DrägerService.
"Accuracy check error"	Accuracy check not successful.	Repeat accuracy check, if fault reappears, call DrägerService.
↑ ↑ ↑	BAC measured exceeds measuring limits.	Apply smaller concentrations.
"Internal System ERROR" with service code	Unit fault.	Remove batteries (2 sec.) and insert them again, if the fault reappears call DrägerService.

Technical Data

Measuring principle

Electrochemical Sensor

Measuring range

Breath alcohol concentration	0.00 to 2.5 mg/L 0 to 2500 µg/L 0.000 to 0.500 g/210L 0 to 250 µg/100mL 0 to 250 µg% (Ethanol mass per breath volume at 34 °C and 1013 hPa.)
Blood alcohol concentration	0.00 to 5.00 ‰ 0.000 to 0.500 % 0.00 to 5.00 g/L 0 to 500 mg/100 mL 0 to 500 mg% 0.00 to 5.00 mg/mL (Ethanol mass per blood volume re- sp. mass at 20 °C and 1013 hPa. By giving this data a conversion fac- tor of 2100 between breath and blood alcohol concentrations is used.)

Display range

The range of indication and the measuring range are identical.
If the result is exceeding the measuring range the display will
show three arrows "↑ ↑↑" as overflow indication.

Ambient conditions

For operation	−5 to 50 °C (23 to 122 °F) 10 to 100 % r.h. (non condensing) 600 to 1400 hPa
For storage	−40 to 65 °C (−40 to 149 °F)

Sample conditions

minimal blowing time	selectable 2 s default
minimal volume	selectable 1.2 L default

Measurement accuracy

Reproducibility with an ethanol standard

up to 0.50 mg/L	0.008 mg/L
≥ 0.50 mg/L	1.7 % of measurement volume the higher volume is valid
up to 1.00 ‰	0.017 ‰
≥ 1.00 ‰	1.7 % of measurement volume the higher volume is valid

Values for other measuring units can be calculated according to
table on page 45.

Sensitivity drift

typically 0.6 % of meas. value/month

Calibration**Intervals before readiness for operation**

1st measurement about 6 s

after terminating a result of a previous test in range:
(breath alcohol unit)

0.25 mg/L	20 s
0.50 mg/L	40 s
1.00 mg/L	80 s
> 1.5 mg/L	120 s

after terminating a result of a previous test in range:
(blood alcohol unit)

0.50 ‰	20 s
1.00 ‰	40 s
2.00 ‰	80 s
> 3.00 ‰	120 s

Intervals for other measuring units selected can be calculated accordingly (see page 45).

Calibration intervaldepends on required accuracy,
typically: 6 months

Permissible calibration range for unit of measuring value:

mg/L	0.160 to 0.750 mg/L
µg/L	160 to 750 µg/L
g/210L	0.035 to 0.150 g/210L
µg/100 mL	16 to 75 µg/100 mL
µg%	16 to 75 µg%
‰	0.35 to 1.50 ‰
%	0.035 to 0.150 %
g/L	0.35 to 1.50 g/L
mg/100 mL	35 to 150 mg/100 mL
mg%	35 to 150 mg%
mg/mL	0.35 to 1.50 mg/mL

Standard calibration value when using Dräger Ethanol-Standard
solution for unit of measuring value:

mg/L	0.476 mg/L
µg/L	476 µg/L
g/210L	0.100 g/210L
µg/100 mL	47.6 µg/100 mL
µg%	47.6 µg%
‰	1.00 ‰
%	0.100 %
g/L	100 g/L
mg/100 mL	100 mg/100 mL
mg%	100 mg%

The use of these standard calibration values corresponds to a conversion factor of 2100 between breath and blood alcohol concentration-values.

Dimensions 140 mm x 80 mm x 35 mm (W x H x D)

Weight approx. 195 g

Electrical supply

Batteries 2 x 1.5 V Type Alkaline
(Mignon, LR6, AA)
approx. 1500 measurements

CE marks Electromagnetic compatibility
(Directive 89/336/EEC)

Order List

Designation and description	Part No.
Dräger Alcotest® 6510 (instrument, 3 mouthpieces, 2 batteries, hand strap and case)	83 17 900
Accessories	
Mouthpieces hygienically tested	
100 pieces	68 10 690
250 pieces	68 10 825
1000 pieces	68 10 830
Alkaline battery (1 piece)	13 35 804
Rechargeable battery NiMH (1 piece)	83 179 12
Case for Dräger Alcotest® 6510	83 179 20
Leather case for Dräger Alcotest® 6510 (belt fixation)	83 17 911
Leather case for Dräger Alcotest® 6510 (yellow) (Instrument can be operated while staying in the pouch)	83 17 931
Calibration accessories	
Mark II Simulator	45 40 088
Ethanol, standard solution 1.0 ‰ (1 x 500 mL)	67 28 838
Calibration set	83 17 927
Calibration gas, ethanol-nitrogen	on request

Dräger Safety AG & Co. KGaA

Revalstrasse 1

D-23560 Lübeck

Deutschland

Tel. +49 451 882-0

Fax +49 451 882-20 80

www.draeger.com

90 23 757- GA 4754.500 de/en

© Dräger Safety AG & Co. KGaA

Edition 04 - 12/2008

(Edition 01 - 04/2004)

Subject to alteration